extróma®

SLIM

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SPARE PARTS CATALOGUE

GENERAL SAFETY RULES

These safety rules are an integral part of the product. Read the information in this manual carefully since it provides important instructions for safe use and maintenance of the system. Keep these instructions in a safe place and ensure that anyone operating the machinery is familiar with them. This product should be used only for the specific purpose for which it is designed: any other use is improper and hazardous.

The manufacturer will not be held liable for damages caused by improper, incorrect or unreasonable use. Do not allow children to play or stand about in the area of action of the machine; do not allow unaccompanied children to use the elevator.

In case of failure or malfunction of the product, disconnect the power switch and do not attempt to repair the machine yourself; contact authorised professional technicians for this purpose. All maintenance and repairs must be done only by professional technicians, authorised for the purpose. To ensure the efficient and correct operation of the plant, observe the manufacturer's instruction regarding scheduled maintenance by authorised technicians; in particular, all safety equipment must be regularly checked. All installation, maintenance and repair work must be registered and the registers made available to the user.

Failure to comply with the above may generate hazards

Contents of the manual

- 1/ Conformity
- 2/ Characteristics and description of the machine
- 3/ Identification plate data
- 4/ Technical service
- 5/ Commissioning
- 6/ Proper and improper use
- 7/ Correct use of the elevator
- 8/ Safety systems
- 9/ Emergency operation by the user
- 10/ Vibration and noise
- 11/ Wiring diagrams
- 12/ Maintenance and inspections
- 13/ Disposal of substances and waste materials

1) CONFORMITY

With the aim of ensuring the highest levels of safety for the user, the design of the machinery and the installation of the SLIM inclined platform have been executed in accordance with the following safety regulations and legislation.

Machine:	European Machinery Directive 98/37 European Electromagnetic Compatibility Directive 86/336 Italian Ministerial Decree 89/236 Wheelchair Platform Lift Standard UNI 9801 Wheelchair Platform Lift Standard TUV 103 - B Machinery Safety Standard EN 292 - 1 - 2 Machinery Safety Standard EN 418 Standard Governing Electrical Equipment of Machinery EN 60204
----------	--

2) CHARACTERISTICS

Components of the plant:

Fixed parts:	rails rail mounts drive unit
Moving parts:	machine body
Control parts:	electrical cabinet control panel

Specifications:

Direction of travel	up/down – electrically driven Speed 0.1m/sec. Capacity 1 person on wheelchair Load 225 daN max (see data plate) Standard duty 30 cycles/hour
Ambient conditions	-10°C to +60°C - max. humidity 70% On-board controls travel direction buttons locking emergency stop button
Floor controls	call button key-switch

The SLIM wheelchair lift is designed for transporting disabled persons on wheelchairs.

The machine travels on two metal rails designed for this application, supported by feet mounted to the stairs and the wall, if present.

The drive unit with electric motor is located at the top end of the rails.

The machine is driven by a metal rope running inside the rail, which is anchored to the upper carriage.

An automatic levelling device keeps the platform horizontal no matter the angle of the rail.

The machine is composed of a fold-away platform with non-slip surface equipped with automatic ramps for connection to the floors, the machine body equipped with control buttons and fold-away safety bars.

The low speed of travel, hold-down controls, collision safety equipment, mechanical safety brake and overtravel terminal enable the unit to be installed on any kind of route, whether straight, curved, etc. in total safety. The floors are equipped with lift call buttons which are key-switch operated (optional).

3) IDENTIFICATION PLATE DATA

- manufacturer	Extrema srl via dell' Industria 2 – 46031 Bagnolo S.Vito (MN) tel. 039 / 0376 / 252443fax 039 / 0376 / 251091
- model	see data plate
- year of manufacture	see data plate
- serial number	see data plate
- load	225 daN max(see data plate)
- capacity	1 person on wheelchair
- speed	0.1m/sec.
-operating voltage	230 V AC 50 Hz
- control voltage	24 V DC
- consumption	2 kW max

4) TECHNICAL SERVICE

Stamp of retailer or authorised service agent

5) COMMISSIONING

After installation, final testing and hand-over to the customer, done by an authorised technician, the vertical platform is put into service by moving the mains switch on the electrical cabinet to ON.
In case of lengthy inactivity, we recommend switching off the plant.
Move the switch to OFF to switch it off.

6) PROPER AND IMPROPER USE

6.1) Proper use

INSTALLATION

Internal or external on existing stairs or route.

The customer is provided with the installation drawings including the mounting points for verification of the strength of the masonry and for building the shafts (if applicable).

WARNING:

Failure to observe the specifications of the installation design drawings can be the cause of hazard.

USE

Transport of only one person on a wheelchair with large rear wheels and small caster wheels at the front, in the direction of travel.

Alternatively, the unit can be equipped with a fold-away seat mounted to the machine body for only one person (optional).

The unit must be operated by persons who are physically and psychologically equipped to do so, and who are familiar with the operation, and the user and maintenance instructions for the device.

WARNING:

If the plant is to be used by a person who is not self-sufficient, it must be operated by a attendant.

6.2) Improper use

INSTALLATION

Do not install the plant in areas subject to the risk of flooding or explosion.

WARNING:

In case of failure to notify the possibility of such events, the manufacturer declines all liability for damages and personal injury.

USE

Use of the plant by unauthorised persons.

Transporting two or more persons on the wheelchair.

Transporting one or more persons upstanding on the platform.

Transporting loads on the platform.

Overloading the machine.

Obstacles inserted between the rail components or placed on the same rail.

WARNING:

Do not put body's parties or objects in the spaces between fixed components ad the ones in motion.

Do not put hands on the rail during the motion of the platform.

WARNING:

**Do not insert the fingers in the upper rail opening:
Danger for squashing!!!**

WARNING:

Do not pour liquids or insert foreign bodies into the holes or slots and do not make unauthorised modifications.

7) CORRECT USE OF THE WHEELCHAIR LIFT

WARNING:

**Use the plant in accordance with the proper use specifications of paragraph 6.
Read this manual thoroughly before operating the system.
Keep this manual in the vicinity of the machine.**

FLOOR CONTROLS

Button control panel located next to the floor halt:

- hold-down travel button
when pressed, automatically moves the machine to the floor in question from any other halt
- key-switch for enabling the controls
turn the key clockwise to enable the call button

BASIC version (manual safety bar operation)

- To close the machine: with the machine halted at the floor and no-one on board, raise the platform to release the lock and rotate the bars downwards; rotate the bars until they are vertically at the bottom with the platform at the top.
In this position the floor controls are enabled
- Enabling the floor controls: turn the key clockwise and hold the call button pressed until the machine arrives, at which point it will halt automatically.
- To open the machine: with the machine halted at the floor, lower the platform; rotate the bars upwards into the open position; the ramp will automatically open towards the floor to allow the passenger to embark.
WARNING: when operated as above, the bar and ramp on the disembarkation side will lock in the safety position to prevent the passenger falling down the stairs. In this position the floor controls are disabled

WARNING:

Check that in the range of operation of the platform and bars there are no obstacles or impediments which may cause damage or personal injury.

AUTOMATIC version (motorised bars and platform)

- To close the machine: with the machine halted at the floor and no-one on board, enable the floor controls by turning the key clockwise; press the button with
- 

the symbol; the platform and bars will move sequentially to the standby position with the bars below and the platform at the top.
In this position the floor controls are enabled

- Moving the machine: hold down the travel button until the machine arrives at the floor, at which point it will halt automatically.

To open the machine:

with the machine halted at the floor, enable the floor controls by turning

the key clockwise; press the button with the  symbol; the platform and bars will move sequentially to the open position; the ramp will automatically open towards the floor to allow access to the platform.

WARNING: when operated as above, the bar and ramp on the disembarkation side will lock in the safety position to prevent the passenger falling down the stairs. In this position the floor controls are disabled

Machine with ramp front access to the platform (optional):

on the AUTOMATIC version only, the motorised front access ramp moves in the same way as the disembarkation bar
bar vertically up = front access open
bar horizontal = front access ramp functions as foot guard
bar vertically down = front access folded away

WARNING:

Check that in the range of operation of the platform and bars there are no obstacles or impediments which may cause damage or personal injury.

ON-BOARD CONTROLS

Controls located on the upper section of the machine body (control board):

- hold-down travel button
when pressed, automatically moves the machine in the indicated direction from the starting floor to any other halt (floor 0 is always the lowest floor)
- red emergency STOP button
when pressed, locks down and instantly halts the machine in position, whatever the direction of travel; to reset the controls, turn clockwise until it releases

Hold down the button until the machine arrives at the floor where it halts automatically and releases the access side bar to allow disembarkation.

For intermediate halts and AUTOMATIC version (motorised bars and platform)

systems with intermediate halts are equipped with a control panel for the passenger's attendant with 4 buttons (up/down; open/close bars).

In this version the bars are opened at the floor manually by pressing the button on the attendant's control panel.

When the machine stops at an intermediate halt, a light or acoustic signal is activated. When the travel button is released, the machine stops and pressing the open button the bar facing the floor raises to enable disembarkation.

WARNING:

**The position of the wheelchair occupant must allow for access to the on-board controls and must be as far away as possible from the unprotected area.
Prevent accidental movement of the wheelchair by applying its brakes.**

WARNING:

Do not introduce parts of the body or other objects into the gaps between the fixed and moving parts.

WARNING:

Do not operate the plant unless there is a person nearby who is familiar with the manual emergency controls.

IN GENERAL

All motorised controls are delayed by 3 seconds and opening the bars vertically upwards disables the drive motor.

To prevent hazardous situations, the operator must be sure that the shaft and stairs to which the platform is installed are adequately illuminated.

WARNING:

The plant may only be used by persons who are authorised, capable of operating motorised machinery and familiar with the use of the plant.

8) SAFETY SYSTEMS

MECHANICAL SAFETIES

Structural dimensions =	designed as required by technical standards
Drive unit =	VSF reducer unit and gears for optimal load distribution
Speed limiter/ mechanical safety brake =	integrated into the mechanical safety brake mechanism, trips in case of excessive downwards travel speed
Mechanical overtravel stop =	stops the upwards travel of the platform
Bars unlockable from outside =	<u>BASIC versions</u> : press the push rods behind the machine's body <u>AUTOMATIC version</u> : engage the manual handwheel in the holes above the machine body

WARNING:

**Falling hazard!
This operation is only for trained persons, who are familiar with rescue procedures.**

ELECTRICAL SAFETIES

Component ratings =	designed as required by technical standards
Safety/control circuits =	low operating voltage
Overtravel cam =	stops the upwards travel of the machine beyond the top floor and disables the system

WARNING:

**For rescue of the passengers, refer to paragraph 9.
To re-activate the plant, we recommend contacting your authorised service centre.**

Cam = trips to disable the plant if the downward speed speed limiter exceeds a set value and mechanical safety brake

WARNING:

To re-activate the plant, you must contact your authorised service centre.

Bar lock control = enables the machine to operate when the bars are in the safety (horizontal) position

Emergency stop button = when pressed, locks down and instantly halts the machine in position whatever the direction of travel; to reset the controls, turn clockwise until it releases

Collision sensors = when pressed, stop the machine travelling; the machine can only start again if the sensors are disengaged. They are located along the sides of the machine body under the platform near to the ramps.

9) EMERGENCY OPERATION BY THE USER

OPERATION FROM THE FLOORS

If, when the call button is released, the platform continues to travel, disable the control by turning the key-switch key counterclockwise to OFF.

WARNING:

If this is not done immediately, the machine will continue to travel upwards/downwards until the safety sensors trip and disable the plant which requires the intervention of an authorised technician to return to normal service. To rescue the occupant, refer to the following paragraphs.

This operation is only necessary when the control circuit is malfunctioning: contact your authorised service centre to have the machine checked over.

ON-BOARD OPERATION

If, when the travel button is released, the platform continues to travel, disable the control by pressing the emergency STOP button.

WARNING:

If this is not done immediately, the machine will continue to travel upwards/downwards until the safety sensors trip and disable the plant. To rescue the occupant, refer to the following paragraphs. This operation is only necessary when the control circuit is malfunctioning: contact your authorised service centre to have the machine checked over.

TO RESCUE THE OCCUPANT

In case of electrical power failure, the travel of the platform can be continued by fitting the manual handwheel into the hole at the base of the drive unit so that it engages the shaft of the electric motor.

WARNING:

Do this only after having switched off the system by turning to OFF the mains switch on the electrical cabinet.

If the platform is blocked along its route, for example due to a power failure, the occupant can be rescued by releasing and opening the bar facing the upper section of the route.

For releasing the bar refer to paragraph 8- "Bars unlockable from outside".

When the bars are unlocked, the platform is disabled; to return to normal operation, move the bar back to the horizontal position.

WARNING: **Falling hazard! Only unlock the bar facing the upper part of the platform's route.**

FREEING UP THE ROUTE (STAIRWAY)

After the occupant has been rescued, the route can be freed up as follows:

BASIC version (manual safety bar operation)

Raise the platform to release the lock and rotate the bars downwards; rotate the bars until they are down with the platform up.

In this position the machine is in service and the floor controls are enabled

AUTOMATIC version (motorised bars and platform)

Fit the manual handwheel into the holes above the machine body and turn it to bring the bars down, then raise the platform.

In this position the machine is NOT in service and the floor controls are disabled

To return the machine to normal operation, move the bars and platform back to the horizontal position.

10) VIBRATION AND NOISE

The body and limbs of the user are subject to low frequency vibrations with very limited accelerations and for short cycles:

These factors are negligible as regards safety.

Noise level measured in the air: **less than 70 dB(A).**

11) WIRING DIAGRAMS

The wiring diagrams are provided as attachments to the documentation.

12) MAINTENANCE AND INSPECTIONS

WARNING: **To ensure an adequate level of safety, observe the specified maintenance intervals and use original spare parts.**

Maintenance is divided into two types:

USER MAINTENANCE

-- Normal cleaning (use a cloth with biodegradable detergent)

WARNING: **Before cleaning the system, switch the machine off by turning the mains switch on the electrical cabinet to OFF.**

MONTHLY INSPECTIONS:

- Check the operation of the control panel enabling keys turning the key to OFF should disable the controls
- Check the operation of the emergency STOP button pressing the button should disable the on-board and floor controls

WARNING: If the machine does not pass these checks disable the plant by turning the mains switch to OFF on the electrical cabinet and immediately contact your authorised service centre.

MAINTENANCE TO BE DONE BY AUTHORISED TECHNICIANS

SIX-MONTHLY INSPECTIONS:


- Check the operation of the mechanical safeties
(wear of the transport rope, operation of the speed limiter/safety brake, mechanical overtravel limit, unlocking the bars from outside, stability of rails and fastenings, grease rails and all rotation axis point).
- Check the operation of the electrical safeties
(overtravel cams, speed limiter and safety brake, sliding contacts, detection of the platform at the floor halt, locking of the bars, emergency STOP, on-board and floor controls)

WARNING: Switch the machine off at the mains switch before working on the electrical cabinet.


13) DISPOSAL OF SUBSTANCES AND WASTE MATERIALS

- The plant does not contain toxic substances in need of special disposal.
- All spare parts, such as cables, cams, etc. in rubber and plastic, should be delivered to authorised collection and disposal centres as provided by established legislation
- Exhausted oils and greases should be delivered to authorised collection and disposal centres as provided by established legislation.

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



extrisma
CODICE
63G300006/c

Instructions for thr insulation tests

- A) Stop the platform between two floors
- B) Switch the main switch OFF
- C) Disconnect the feeding from inverter and take out the fuse FU4 body
- D) Disconnect the ground cable (earth) of the feeding cable
Operation circuit: check the insulation to earth of "90" terminal
- E) Manually switch the electromagnetic power switches KM1 and KM2 on
Traction motor circuit: check the insulation to earth of "U1", "V1", and "W1" terminal
Finally check the reciprocal insulation between terminals "90" and "U1", "90" and "V1", "90" and "W1"

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SIMB.	DESCRIPTION	RAMO
AP1	MAIN ELECTRONIC PLC PRINT	20-->61
AP2	CONTROL ONBOARD PRINT	68-->99
B1	EMERGENCY SIREN	55-->59
C1	CAPACITOR 4700µF 50V	11
E1	4 CHANNELS RECEIVER	26-->30
E3	FLOOR CONTROL	40-->45
E4	ON-BOARD CONTROL	76-->77
E5	CONTROL FOR ATTENDANT	70-->75
FU1(A,B)	MAIN LINE FUSES 2X10A	03
FU2	TRANSFORMER FUSE 1A	05
FU3	SUB-CIRCUIT FUSE 10A	13
FU4	AUXILIARY FUSE 24 VDC 10A	13
KM1	MAIN CONTACTOR 1	00-03-22-23
KM2	MAIN CONTACTOR 2	00-03-24-35
M1	ELECTRIC MOTION ENGINE	03
M2	ELECTRIC ENGINE FOR RIGHT ARM	89
M3	ELECTRIC ENGINE FOR LEFT ARM	81
M4	ELECTRIC ENGINE PLATFORM MOTION	91
M5	ELECTRIC ENGINE SIDE ACCESS RAMP	95
QS1	MAIN SWITCH	03
SA1	FLOOR CONTROL KEY SWITCH	41
SA2	RIGHT/LEFT SELECTION PLAN	83-87
SBE	STOP-PUSHBUTTON ON BOARD	65
SB0	FLOOR UPWARD PUSHBUTTON	41
SB1	FLOOR DOWNWARD PUSHBUTTON	42
SB2	FLOOR PUSHBUTTON CLOSING	43
SB3	FLOOR PUSHBUTTON OPENING	44
SB90	ON-BOARD UPWARD PUSHBUTTON	76
SB90'	FLEXIBLE ON-BOARD UPWARD PUSHBUTTON	71
SB91	ON-BOARD DOWNWARD PUSHBUTTON	77
SB91'	FLEXIBLE ON-BOARD DOWNWARD PUSHBUTTON	72
SB92	ON-BOARD PUSHBUTTON "ARMS OPENING"	74
SB93	ON-BOARD PUSHBUTTON "ARMS CLOSING"	75
SQ1	CLOSED RIGHT ARM MICROSWITCH	84
SQ2	OPENED RIGHT ARM MICROSWITCH	85
SQ3	HORIZONTAL RIGHT ARM MICROSWITCH	78-83-86
SQ4	CLOSED LEFT ARM MICROSWITCH	87
SQ5	OPENED LEFT ARM MICROSWITCH	88
SQ6	HORIZONTAL LEFT ARM MICROSWITCH	78-88-89
SQ9	25° PLATFORM SETTING MICROSWITCH	90
SQ10	OPENED PLATFORM MICROSWITCH	92
SQ11	CLOSED PLATFORM MICROSWITCH	93
SQ12	FOLDED PLATFORM CHECK MICROSWITCH	78-79
SQ13a..h	LEFT/RIGHT ANTISHOCKING MICROSWITCH	65-67
SQ14a..d	PLATFORM BOTTOM ANTISHOCKING MICROSWITCH	67
SQ15	SIDE ACCESS SAFETY SETTIG MICROSWITCH	78-94

SIMB.	DESCRIPTION	RAMO
SQ16	OPENED SIDE RAMP MICROSWITCH	95
SQ17	CLOSED SIDE RAMP MICROSWITCH	96
SQ18	30° PLATFORM SETTING MICROSWITCH	85
SQ19	Parachute switch	25
SQ29a-b	Safety limit switch	25
SQ39	SLOW DOWNING MICROSWITCH	69
TC1	Auxiliary transformer	07
VC1	BRAKE ENGINE DIODE BRIDGE RECTIFIER	10
VC2	DIODE RECTIFIER FOR AUXILIARY	10
VD	Diode	23-36
X1	ELECTRIC BOX MAIN TERMINAL BOARD	--
X2	INSTALLATION DEVICE CONNECTION	--
XPS	CRAWLING CONTACTS FLOOR STOP	82-83-84-87
YB1	MOTION ENGINE ELECTROBRAKE	02
U1	INVERTER FOR TRANSLATE MOTORS	03-25-->31

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


SIMB.	DESCRIPTION	RAMO
SQ19	PARACHUTE MICROSWITCH	25
SQ29a-b	SAFETY LIMIT MICROSWITCH	25
SQ39	SLOW DOWNING MICROSWITCH	49
TC1	AUXILIARY TRANSFORMER	07
VC1	BRAKE ENGINE DIODE BRIDGE RECTIFIER	01
VC2	DIODE RECTIFIER FOR AUXILIARY	09
VD	DIODE	23-36
X1	ELECTRIC BOX MAIN TERMINAL BOARD	--
X2	INSTALLATION DEVICE CONNECTION	--
XPS	CRAWLING CONTACTS FLOOR STOP	62-63
YB1	MOTION ENGINE ELECTROBRAKE	02
U1	INVERTER FOR TRANSLATE MOTORS	03-25-->31

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VIETATE LE RIPRODUZIONI NON AUTORIZZATE REPRODUCTION NOT PERMITTED ALL RIGHTS RESERVED